





Protocol for setting up SPA data collection using BIS method.

Before starting make sure that microscope is working correctly. On the microscope computer TEM Imaging Analysis (TIA), Microscope software launcher, FEI-SEM server, and Gatan remote TEM are running. On the Gatan K3 computer make sure digital micrograph and SerialEM are running. Check network connections on the Gatan K3 computer.

This document is for core staff only and researcher that have been trained to use the microscope.

Typical work flow for a single-particle data collection using beam-image shift






1. On the Gatan K3 computer make a new folder in the X: directory to save your data. **Make sure that the X: drive has enough free space if it is full delete old data.** Make sure the data has been fully transferred to the longleaf cluster or the NAS before deleting the data.

Windows File Manager

“This PC” → Capture Data (X) → K3Data → [new session file: date_initials]

ie 20250502_JDS_CBL



At this point you are ready to start the session.

2. Start the Autoloader Inventory on the microscope computer.

If the cassette slot is occupied by a cartridge make it blue using the mouse clicker, if empty make dark grey.

Check temperature of the autoloader and column before starting the Inventory.

Check “Turbo Auto Off” on the Autoloader Option tab (blue arrow)

The screenshot shows the Autoloader (User) interface with several panels and callouts:

- Autoloader (User) Panel:**
 - Cassette:** Dock (yellow), Undock (grey).
 - Cartridge:** Slots 1-12. Slots 1, 2, and 11 are blue. Slot 12 is dark grey. Slots 3-10 are light grey.
 - Buttons:** Edit Slot State (yellow), Stage (grey), Load (grey), Unload (yellow).
 - Status:** (empty).
- Options Panel:**
 - Vacuum:** Turbo Auto Off (default) (selected), Turbo Always On (radio button).
 - Control:** Initialize (yellow), Loader-cycle (grey), Vacuum On (yellow).
 - Inventory:** (yellow button, highlighted by a purple arrow and callout).
 - Cassette Unlock:** (grey).
 - Leave cartridge on CompuStage:** (checked checkbox).
- Temperature Control Panel:**
 - Status:** All Nitrogen Temperature.
 - Dewar levels:**


Autoloader	39 %	3 h 41 min
Column	63 %	8 h 33 min
 - Temperatures:**

Docker	121.2 K	-151.9 °C
Holder	79.4 K	-193.8 °C
Cassette gripper	120.4 K	-152.7 °C
Cartridge gripper	125.4 K	-147.8 °C
Autoloader Dewar	78.8 K	-194.4 °C
Column Dewar	78.5 K	-194.6 °C
- Vacuum (Supervisor) Panel:**
 - Status:** All Vacuum (Closed) (red text).
 - Accelerator:** 1 Log
 - Column:** 1 Log
 - Autoloader:** 27 Log
 - Detection Unit:** 16 Log
 - Buffer tank:** 56 Log
 - Backing line:** 65 Log
 - Buttons:** Col. Valves Closed (yellow), Empty Buffer (grey).
- State Panel:**
 - Temperature State:** Go to: (dropdown), At: 6/18/2025 1:08:25 PM, Set: None.
 - After Conditioning:** Go to: All Nitrogen Temperature (dropdown), At: 6/19/2025 1:08:25 AM, Set: None.
- Settings Panel:**
 - Pressure unit:** Log (dropdown).
 - Pressure color highlighting:** Default (dropdown).
 - Detect. Unit Vented:** (grey button).

Callouts:

- 2nd Left Mouse Click to change color (slot state):** Points to slot 12.
- 1st Click Edit Slot State:** Points to the Edit Slot State button.
- Click inventory to start:** Points to the Inventory button.
- Wait for autoloader to cool to -155 (should be green not red) before inventory:** Points to the temperature table.
- Blue arrow:** Points to the Turbo Auto Off option.

Click Inventory to start (should take ~15 mins)



At this point the Autoloader inventory is complete with no errors on the microscope, as shown in this image.

The screenshot displays the CryoEM Core software interface, divided into four main sections:

- Autoloader (User):** This section contains controls for the autoloader. It includes a "Cassette" section with a "Dock" button and an "Undock" button. Below this is a "Cartridge" section with a list of slots (1-12) and their corresponding labels: CBL1, CBL2, and Lacey. Each slot has a green checkmark icon. At the bottom of this section are buttons for "Edit Slot State", "Stage", "Load", and "Unload".
- Options:** This section contains settings for the autoloader. It includes a "Vacuum" section with radio buttons for "Turbo Auto Off (default)" and "Turbo Always On". Below this is a "Control" section with buttons for "Initialize", "Loadercycle", and "Vacuum On". There is also an "Inventory" button. At the bottom, there is a "Cassette Undock" section with a checked checkbox for "Leave cartridge on CompuStage".
- Temperature Control:** This section displays the status of the temperature control system. It includes a "Status" section with "All Nitrogen Temperature" and "Dewar levels" (Autoloader: 37 %, 3 h 29 min; Column: 62 %, 8 h 27 min). Below this is a "Temperatures" table showing the current and setpoint temperatures for various components.
- State / Filling:** This section displays the current state of the system. It includes a "Temperature State" section with a "Go to" dropdown menu and a date/time selector (6/18/2025 1:08:25 PM). Below this is an "After Conditioning" section with a "Go to" dropdown menu and a date/time selector (6/19/2025 1:08:25 AM).

Temperatures		
Docker	108.6 K	-164.5 °C
Holder	79.4 K	-193.8 °C
Cassette gripper	95.7 K	-177.5 °C
Cartridge gripper	111.1 K	-162.0 °C
Autoloader Dewar	78.8 K	-194.4 °C
Column Dewar	78.5 K	-194.6 °C

Optional do gain reference during the Inventory and direct alignment on evaporated aluminum or another TEM grid (not UltraAufoil).

3. On SerialEM open the current Setting file for SPA.

“Settings” → 20250325_SPA_45k” (Yellow arrow)


Or go to C:\Program Files\SerialEM\

Load one of the setting files

20250325_SPA_45K (normally use this – pixel size 0.8737 Å)

20250427_SPA_73K (Pixel size 0.54925 Å)

20250427_SPA_22K (Pixel size 1.7738 Å)




Make sure Low-Dose Mode is checked in the Low Dose Control Panel (Blue arrow).

Check Record, View, Focus, and Search make sure the beam centered and camera settings are correct for each.

To do this first in the microscope alignment tab select Beam shift and use the multi-function XY knobs on the control panels to center the beam in Record mode. Then select “Additional beam shift Set” in the Low-Dose Control panel to move the beam in View, Focus and Search” using the multi-function XY knobs (do not select beam shift).

To change the beam setting or magnification ect...Click “continuous update” on the Low-dose Control Panel. Normally you will not need to do this.



4. Run the script DoseCal to write to the SerialEM log file what the imaging conditions are.


“Script” → “Edit” → “Script1: DoseCal”

Navigate to a broken grid square or unload grid (preferred). Go to microscope computer and under the Autoloader tab to unload the cartridge.

Go to record in low dose make sure the 70 micron condenser aperture is inserted and centered. Make sure the imaging and camera parameters are exactly what you are using to collect data.

Run “DoseCal” → pop-up: navigate to **[new session file: date_initials]** → “Ok”


{ It writes CS, Voltate, flux, total dose and pixel size to the log file}



If needed using a Cross-replica, evaporated aluminum, or carbon coated Quantifoil grid align the microscope.

5. Montage each of the grids in cassette using script Map2Grids.

Before running script read the header and follow the instructions



Read the script header for instructions.

```
MacroName Map2Grids
# Script map2grids modified from Wim Hagen
# Open Navigator, Setup Full Grid Montage 6 X 10 overlap 23% for UF 15% for QF
# Low-Dose Search should be 62X Mag
# Have the C2 apt 150um inserted
# Complete inventory
# Modified by JDS 20200715
```

“Navigator” → “Open” to open the Navigator

“Script” → “Edit” → “2” (“Script2: Map2Grids”)

Edit the script Parameters

Edit “SetDirectory X:\k3data\” [new session file: date_initials]

Edit “cat”: each cartridge you want to montage

Edit “name” name of file or sample


```
##### Parameters #####
# set directory for saving images
SetDirectory X:\K3Data\20250325_JDS_CBL

cat = { 11 10 9 8 7 6 5}      # change this which cartridge you like to montage
name = { apo1 apo2 apo2 apo4 apo5 apo6 apo7 apo8} # change this to the name of the sample
##### Parameters #####
```

Go to “Search” in low dose and make sure it is set to 62X with the C2 150 micron condenser aperture inserted.


If Search is a different magnification, change it by checking continuous update in the Low Dose Control and set the magnification to 62X. Alternative select the 62X 4uP in the imaging state panel.

Check the Search beam on the FlueCam on the microscope. *If needed center the beam using by checking the “additional beam shift” in low-dose panel and beam-shift using the multi-function XY knobs. If needed center the aperture in the Tab Aperture, click Adjust, and center using the multifunctional ZY knobs.*



In SerialEM

“Navigator” → “Montaging and Grids” → “Setup Full Montage”



In the pop-up click “Ok” normally you don’t need to change anything.

If needed, adjust “number of pieces” (default 6 x 10, full stage limit)

Check camera(K3), Mag (62), Bin (1), Overlap (20%)

“Save As” → Navigator [new session file: date_initials] → [fullgridmontage] → “Save”

“Script2: Map2Grids” → Run (takes about 15 minutes per grid).

Alternative you can montage a single grid on the stage, as described above but don’t run the script just Click “Start” on Montage Control Panel.


It will take ~ 15 mins to montage one grid.

At this point the montaging is done, and the maps are saved in the navigator.


6. Collect grid square montages at each grid square you want to collect data at.

On SerialEM zoom into Full Grid Montage and add point to the center of the grid squares

Navigator → Add Points → stop adding points → mark points Acquire (A)



Find eucentric height (once at one grid square) using stage Wobbler on the microscope computer and “Update Z” for Acquire points (Purple arrow).



**Change Search to 210 X magnification.**


Imaging state → select “Square Search SPA” (mag 210)

Or check continuous update on low-dose panel and change magnification of Search. Uncheck continuous update when done.


Note for 200 mesh grids you will need to use mag of 115 not 210.

Check camera parameters in Search.

Camera → Parameters → “square” → “OK”




Take an image of the grid square and copy the image to Buffer M and N in the Buffer Controls pannel (Yellow Arrow).




Open new montage file

File → Montage Set up

Check camera(K3), Mag (210), Bin (1), Overlap (20%), Check move stage instead of shifting image, check make map from each montage if navigator is open, check close file when montage is done, check use Search parameters in low dose mode.





Save As → Navigator (X:\k3data to [new session file: date_initials] → [MMM] → “save”

Navigator → acquire at items → mapping → “GO” (1.5 min per point)

Acquire at Items

Parameters for: **Mapping** **Final Data**

Primary Action

- ☒ Acquire and save image or montage
 - ☒ Make Navigator map
 - ☐ Skip saving to file
- ☐ Acquire multiple Records Setup
- ☐ Acquire tilt series
- ☐ Run script SquareCenter

Primary Action-related

- ☐ Cycle defocus target in autofocus, drift wait & TS
From -2 to -1 um in 10 steps
- ☐ Adjust beam tilt & astig for image shift
- ☐ Do early return with sum of 0 frames
- ☐ Use map hole vectors for shifts (no adjustment)
- Acquire non-montage images in Low Dose with:
 - ☐ Record
 - ☐ View
 - ☒ Search

MONTAGES will be saved into the file:
MMM.mrc

Task-related options

- ☐ Do hybrid Realign to Item/Align to Template
- ☐ Relax stage after final movement
- ☐ Skip Z moves in initial move and Realign
- ☐ Skip stage move to item if possible

General options

- ☐ Do subset from index 1 to 4
- ☐ No message box when errors occur
- ☒ Close column valves at end
- ☐ Retract cameras at end
- ☐ Send email at end
- ☐ Hide unselected options

Tasks before or after Primary (Push button to show controls)

Task	Frequency	Setup
<input checked="" type="checkbox"/> Manage Dewars/Vacuum	Every item	<small>Setup</small>
<input type="checkbox"/> Coma-free Alignment	Every item	<small>Setup</small>
<input type="checkbox"/> Correct Astigmatism	Every item	<small>Setup</small>
<input type="checkbox"/> Condition Phase Plate	Every item	<small>Setup</small>
<input type="checkbox"/> Update Dark Reference	Every item	
<input type="checkbox"/> Rough Eucentricity	Every item	
<input type="checkbox"/> Eucentricity by Focus	Every item	<small>Setup</small>
<input type="checkbox"/> Autocenter Beam	Every item	<small>Setup</small>
<input type="checkbox"/> Realign to Item	Every item	<small>Setup</small>
<input type="checkbox"/> Cook Specimen	Every item	<small>Setup</small>
<input type="checkbox"/> Fine Eucentricity	Every item	
<input type="checkbox"/> Align to Template	Every item	<small>Setup</small>
<input type="checkbox"/> Autofocus	Every item	<small>Setup</small>
<input type="checkbox"/> Wait for Drift	Every item	<small>Setup</small>
<input checked="" type="checkbox"/> Run Script before Action	Every item	

Script to run: SquareCenter

Primary Action (Acquire Image) occurs

- ☐ Run Script after Action Every item, after
Script to run: # 11
- ☐ Hole Finder & Combiner Every item, after Setup

General controls for task: Correct Astigmatism

How often to run task:

- ☒ Every 1 items
- ☐ Only at start of group
- ☐ Only at end of group
- ☐ Every 15 minutes
- ☐ When moved > 40 um

Position in sequence:


Move Up Move Down

☐ Run after Primary Action

☐ Run at nearest item with ☒ label ☐ note starting with

☐ Hide unused tasks Restore Default Order


Save Parameters to File Read in Parameters **GO** Postpone Cancel ?




7. Setup Multishot Record

Take search image on square with visible Quantifoil holes

On the Navigator panel Add 5 points in dice pattern as shown.




Take view of center point. If point is away from center hole (# 8 as shown), place marker in center of hole, navigator → shift to marker → “OK”



Move stage to the center point.


Navigator panel → “Go to XY”



Take a view, check that it's centered (to center left mouse click and drag to move the stage over the center of the hole).

Navigator → montage and grids → Setup Multishot Parameters

Click “Four corners of regular pattern” (Multiple Record Setup - Pink Arrow) to begin



Navigator

Select 1st point (#9 as shown -- Blue) in multishot → IS to XY (Green Arrow).

Take view.

If not centered, right click and drag to center.

Save Image shift (Multiple Record Setup).

Navigator.

Select 2nd point (# 10 as shown) → IS to XY.

Take view.

If not centered, right click and drag to center.

Save Image shift (Multiple Record Setup).

Navigator.

Select 3rd point (#11 as shown) → IS to XY

Take view

If not centered, right click and drag to center

Save Image shift (Multiple Record Setup)

Navigator

Select 4th point (#12) → IS to XY

Take view

If not centered, right click and drag to center

Save Image shift

Take Search and View image to check pattern over holes


Multishot Record → Change to 7 by 7 (or 5 by 5)

8. Set position of Focus beam

In Low Dose Panel, Define position of area focus → on View image left click on carbon film between Quantifoil holes so the yellow circle in over the carbon.

Set “Position on tilt axis” for R1.2/1.3 that is normally -1.8 or 1.8 microns.

Then click None on the Low Dose Panel to switch back.



At this point Focus beam is over the carbon and not illuminating the hole.



9. Check Focus beam so it does not overexpose sample

In Low Dose panel, go to View to move the stage over a hole in SerialEM by taking images or using the microscope control panels. To do this follow these instructions....

On the Microscope

Insert Flue Cam.

Make sure it is centered over the hole if not adjust the stage XY position.

On SerialEM.

In Low Dose Panel, Select Focus, and click additional beam shift check.

On Microscope Flue Screen.


Adjust position of the beam using the XY multifunction knobs so the beam is centered.

On SerialEM.



In Low Dose panel, uncheck additional beam shift.

At this point you have confirmed that Focus beam is over the carbon and not illuminating the hole.

10. Check Offset to View and Search so Record, View and Search are over the same



area.



When aligned the Record or Preview, View and Search will look like this (above)

On the Navigator (or using the FlueCam) add point or marker on an ice crystal go to XY (as shown above).

Go to Record, search for ice, use Flu Cam on microscope computer, move stage to the ice, take a Preview or Record.....if the image are not over the same area set offset for View and Search in the Low-dose pane.

To do this....

Select offsets for View (Green arrow), take View image.

Uncheck move stage for big mouse shifts (Purple arrow).


Right click and drag to correct place on ice → “set” in lose dose control (Blue arrow)

Re-check move stage for big mouse shifts (Purple arrow).


Select offsets for Search, take Search image (Green arrow).

Uncheck move stage for big mouse shifts (Purple arrow).

Right click and drag over ice → “set” in Low Dose panel (Blue arrow).




Re-check move stage for big mouse shifts (Purple arrow).



Take a Search, View and Preview image, the ice should be in the center of each image.

At this point the Search, View and Record are over the same area of the sample



11. Copy image of handsome Quantifoil hole to buffer P


Go to a Quantifoil hole that has no ice contamination.


Take a View Image of it.

Crop image in Camera parameter if needed.

Process → reduce image (factor 3) → “ok”

Copy to buffer P (Buffer Control).





12. Find-holes on each Grid Square Map.

Double click on a Grid Square (MMM) map in navigator panel to open it.

Add points to the map using the hole finder.



Navigator → Montage and Grids → Find holes in regular grid.

Adjust hole size and periodicity as needed (Orange arrow).


To measure holes size Shift + Left mouse click, drag to draw line over holes. look at log file for values.

Click “Find holes” (Green arrow).


Adjust sliders in the Hole Finder so the acquire points turn (pink) and exclude holes (Blue or cyan, as shown).



Click “Make navigator point” (Hole Finder – previous page)




Navigator → montage and Grids → combine points for multishot



Click “Combine Points”

Do this for each MMM Map.

At this point each of the Grid Square Maps (that are good) have Acquire Points over the holes for data collection.



13. Setup data collection

Open and edit data collection script, normally don't need to change anything. The defocus range might need to be adjusted from experiment to experiment.

Normally set

Touch = 0 {to move the beam away from the center of the hole}

Fixobj = 0 { to correct Obj Astigmatism}

TD_low = 0 to - 0.5

TD_high = -0.5 to -1.5

“Script” → “Edit” → “4” (“Script4: DataCollection_V6”)

```
##### SETTINGS #####
maxholeshift      = 200    # max beam/image shift after hole centering in nanometers
driftcrit          = 2.0    # A/sec
driftinterval     = 10     # seconds of delay between drift measurements
drifttimes        = 7     # max no of drift measurements
TD_low            = -0.00   # low defocus
TD_high           = -0.50   #high defocus microns
delta             = 0.25   # defocus step microns 0.1
touch             = 0      #set to 0 for off
bmsft             = 0.300
stagex = 1        # For QF grids set to hole size, if UltraAufoil set to 0
fixobj = 1        # set to 1 to enable and 0 to disable
```

Adjust “Script5” (Multishot Threshold) change the “numMultishot” value. As an example for a 7x7 multishot set to 25, for 5x5 multishot set to 20.

```
MacroName Multishot_Filter_Threshold
# Written JDS on 4/10/22
# Run this script before action when using Data_Collection_V3
# change "numMultishot" before running, will skip points with have less Multishot holes than this
##### Settings #####
numMultishot = 25 # Threshold for acquiring multishot, less than this the item will not be acquired
#####
```

Camera → Parameters → record

As shown in the image.....

Check save frames.


Uncheck align frames.

Set File options → update basename to date → OK


Set folder → navigate to (X:K3Data\ [new session file: date_initials])



Make folder [Data] navigate to [Data] → OK



Navigator → acquire at items → final data → run script “DataCollection_V6”, run script before action Multishot_Filter_Threshold → GO (couple hours)



As shown in this image.



At this point data is being collected, make sure to check the Gatan K3 computer.

14. Run CryoSPARC Live to assess data quality.

Typically, the MaxCTF fit is $< 5 \text{ \AA}$

Defocus range 5000 – 20000 \AA

Total motion < 50 pixels

Astigmatism is 0-600 \AA

Particle are visible and embedded in thin ice.

If the data quality is poor check the microscope, data collection parameters, serialEM setup.

15. If you are collecting data O/N or over weekend Run Argus2.py on workstation #2.

On WKS1 open the command terminal

CD to home directory (/home/cryosparc_user)

```
$ cd ~
```

Now run argus.2 and enter the directory where the data is

```
$ python3 ./argus2.py
```

```
Enter /directory/to/monitor : /mnt/k3data/20250603_JDS/Data
```

```
Do you wish to send additional notifications? Enter (y/n): n
```

```
Skipping additional notifications.
```

```
Argus is now monitoring /mnt/k3data/20250603_JDS/Data
```

If you type Y, You can add your cell phone or email to get a notification once data collection stops.

16. Monitor the data collection for 15-30 mins to make sure there are no errors on the microscope side. If you are collecting data O/N please make sure the data quality is good, argus2 is running, and that there are no errors on the microscope or camera.

The 2D Classes averages should have some ultrastructure, shape or elements of 2nd structure. Depending on the sample wait until 100-200 movies have been collected or ~5,000 particles in the class averages.

Before you leave the core check...

1. LN2 tanks on the microscope.
2. No errors on the microscope, serialEM, or workstation.
3. Walk in the microscope room make sure nothing is out of place or missing.
4. Check the Gatan K3 computer make sure it has enough disk space.