



Clipping Handout

- i. Know what grids you want to clip
- ii. Reserve clipping station in iLab
- iii. Alert Cryo-EM staff to reservation if any assistance, supervision, or grid retrieval will be needed

1) Acquire all equipment, tools, and supplies

- | | |
|---|--------------------------|
| a. Outer plastic piece | i. The 8 clipping tools |
| b. Blue foam piece | j. New rings |
| c. Inner plastic piece | k. New clips |
| d. Silver metal part of clipping station (pedestal) | l. Short, blunt tweezers |
| e. Gold metal part of clipping station (dial) | m. Long, fine tweezers |
| f. Dissection forceps | n. Autoloader tweezers |
| g. Screwdrivers | o. Autoloader boxes |
| h. Ethanol waste | p. Lamp |
| | q. Hairdryer |
| | r. Kim Wipes |

2) Set-up Clipping Station

- a. Ensure that all the equipment parts and pieces are completely dry, use hairdryer and Kim Wipes as needed
- b. Assemble equipment
- c. Cool fully assembled clipping station with liquid nitrogen
- d. Pick, label, and place autoloader grid boxes
- e. Place unclipped grid boxes
- f. Loosen unclipped grid box screws

3) Clip Grids

- a. Use tweezer handle to spin dial to fully reveal groove for ring
- b. Use **BLUNT** tweezers to place ring in the groove
- c. Spin dial to cup ring in a tunnel
- d. Pick up grid with **LONG** tweezers and float grid into the tunnel over the ring
- e. Tap station as needed to center grid in the ring
- f. Use **BLUNT** tweezers to transfer clip to squeeze ring container pad
- g. Use **BLUNT** tweezers to deeply grab the center of clip
- h. Wedge clip fully into clipping tool
- i. Use clipping tool to stamp clip into ring
- j. Use **AUTOLOADER** tweezers to place ringed grid into autoloader box

4) Clean Up

- a. Tighten grid box screws
- b. Fill blue foam dewar with LN2 and quickly transfer grid boxes into it
- c. Place pedestal and dial on the heating block to evaporate
- d. Pour out the clipping station LN2
- e. Place the clipping station in the dehydrator to evaporate, separating all parts
- f. Return all tools to their original locations

5) **FILL OUT CLIPPING STATION NOTEBOOK**