



This training is to ensure that the Vitrobot does not become damaged, samples are properly prepared, and that users remain safe.

Vitrobot Handout

- i. Plan the number of samples, sample concentrations(s), and Vitrobot settings
 - ii. Reserve Vitrobot in iLab
 - iii. Alert Cryo-EM staff to reservation if any assistance or supervision will be needed
 - iv. Bring any supplies and/or tools you wish to use instead of what is already provided by the Core
 - v. (recommended) Bring a facemask to protect samples from your breath
- Set-up Vitrobot
 - Turn on back switch
 - Attach humidifier, syringe in 20-30mL of distilled H₂O, and syringe out the excess
 - Using gloves and flat forceps, remove old blot papers and replace with new – facing the lipped edge towards the blotting pads
 - Set Vitrobot Settings: (recommended)
 - Temperature 22° or 4°
 - Humidity (Off/On/Manual) On
 - Humidity (%) 95%
 - Blot Time BT_i 2, 3, or 4
 - Wait Time WT 0
 - Drain Time DT 0
 - Blot Force BF -10
 - Blot Total BT_o 1
 - Skip Application leave unchecked
 - Use Footpedal check
 - Humidity off During Process check
 - Skip Grid Transfer check
 - Autoraise Ethanelift leave unchecked
 - Fill Metal Dewar with Liquid Nitrogen
 - Assemble Vitrobot Pot:
 - Styrofoam pot, free of grids
 - Styrofoam ring, **PLACE WITH LIP AT THE TOP**
 - Ethane cup
 - Pedestal
 - Spider
 - Cool Vitrobot Pot with Liquid Nitrogen
 - Plasma Clean Grids
 - Add Ethane-propane to Ethane Cup
 - Check that no LN2 is in the ethane cup, if so, soak up with a Kimwipe
 - Place ethane-propane needle nozzle in the ethane cup
 - Open ethane-propane needle valve, **DO NOT REMOVE YOUR HAND UNTIL DONE**
 - Fill ethane cup until touching under the spider, pulling out needle nozzle while topping it off
 - Close ethane-propane needle valve
 - Wait until the surface of the ethane-propane measures quickly to below -190°C
 - **REMOVE SPIDER w/ tongs, don't touch w/ your hands, it is very cold**
 - Label Grid Boxes (recommended: initials, date, alpha-numeric [e.g. CBL 040224 A])

- Dry Run Vitrobot Settings to Test that It's Working Properly (see steps below)
- Operate Vitrobot
 - Tap "Place New Grid"
 - Pick up grid with Vitrobot forceps **AND CLAMP TO THE FIRST MARK**
 - **PUT VITROBOT FORCEPS ON VITROBOT ROD AND CONFIRM THAT IT'S CENTERED**
 - Tap "Next" or footpedal to move arm up
 - Place Vitrobot pot
 - Tap "Next" or footpedal and wait for Vitrobot pot to move into place
 - Tap "Next" or footpedal to lower grid for applying sample
 - Pipet sample onto grid (recommended: 3μL)
 - Tap "Next" or footpedal for Vitrobot to apply WT, BT_i, BT_o, and DT settings and then plunge
 - Spin the Vitrobot pot as needed for ease of grid transfer
 - Confirm that the window is open over the grid box slot
 - Remove the forceps from Vitrobot arm, keeping grid under ethane-propane
 - Quickly dunk the grid into LN2 nearest to open cryo-box slot
 - Quickly place grid into cryo-box slot
 - Use forceps to spin lid so that the window to next grid slot is open
 - Place forceps onto heating blocks to warm up and melt
 - When melted, use Kimwipe to dry forceps or use the hair dryer
- Clean Up
 - Tighten grid box screws
 - Fill foam dewar with LN2 and quickly transfer grid boxes into it
 - Place Vitrobot pot under the snorkel to evaporate
 - **TAP "Exit" ON VITROBOT**
 - Wait for Vitrobot to shut down, and then turn off back switch
 - Detach humidifier, tip upside down into foam box, and syringe out H₂O
 - Return all tools to their original locations
- **FILL OUT VITROBOT NOTEBOOK**

Providing the resources, training, and technical assistance in all aspects of cryoEM since 2019. Written by Clara Lenger