
Optimal Flow Rate and cell concentration for Attune NxT

Sample Flow Rate	Maximum sample concentration for analysis (<u>35,000 evts/sec</u>)	Maximum sample concentration for accurate cell counts (<u>8,000 evts/sec</u>)	Description
1000 $\mu\text{L}/\text{min}$	2.1×10^6 cells/mL	4.8×10^5 cells/mL	Particles > 4 μm
500 $\mu\text{L}/\text{min}$	4.2×10^6 cells/mL	9.6×10^5 cells/mL	
200 $\mu\text{L}/\text{min}$	6.7×10^6 cells/mL	1.5×10^6 cells/mL	Particles > 2 μm
100 $\mu\text{L}/\text{min}$	1.3×10^7 cells/mL	3.0×10^6 cells/mL	
25 $\mu\text{L}/\text{min}$	5.4×10^7 cells/mL	1.2×10^7 cells/mL	Small particles < 2 μm
12.5 $\mu\text{L}/\text{min}$	1.0×10^8 cells/mL	2.4×10^7 cells/mL	Best resolution from background for dimly positive assays

*This information is only to serve as a guide. The biology of your cells and the quality of your data should govern the flow rate and concentration of your experiment