



UNC Flow Cytometry Core Facility

Title:

EMD Millipore Amnis ImageStreamX MkII INSPIRE User Training

Classification:

User Training

Effective Date: 01/03/2017

Revision Date: 01/25/2023

ID: UT SOP004.2

Page 1 of 4

Table of Contents

1. Fluidics.....	2
2. Startup.....	2
3. Acquisition	2
4. Shutdown.....	3
5. Data Analysis Tips.....	4
6. References.....	4
7. Revisions.....	4



UNC Flow Cytometry Core Facility

Title:

EMD Millipore Amnis ImageStreamX MkII INSPIRE User Training

Classification:

User Training

Effective Date: 01/03/2017

Revision Date: 01/25/2023

ID: UT SOP004.2

Page 2 of 4

1. Fluidics



- **SpeedBeads:** monitor and synchronize the flow of the sample and maintain focus and core tracking
- **Sterilizer** = FACSClean
- **Cleanser** = Coulter Clenz
- **Debubbler** = 70% Isopropyl alcohol
- **Rinse** = Deionized Water (diH₂O)
- **Sheath** = 1X PBS (Phosphate-buffered saline), 0.1 nm filtered

2. Startup

- First Person of the day will **Startup**- Initializes fluidics by flushing sheath and loading beads ~14min
- Select "Run **ASSIST** (**A**utomated **S**uite of **S**ystemwide **I**mage**S**tream **T**ests) after initialization" Calibration and testing using SpeedBeads -20min
- If instrument is already started, **Start Fluidics** to ensure Focus and Centering is valid using the SpeedBeads

3. Acquisition

- File > **Load default template** or experiment template.
- **File Acquisition:** set path and number of objects to collect
 - a. Always save to the Desktop initially then copy to the J: Drive ONYEN folder
- **Illumination:** turn on the appropriate lasers.
 - a. Begin with all lasers at max output (hover over input box to view range of power available per laser)
 - b. Keep SSC power between 3-5 mV



UNC Flow Cytometry Core Facility

Title:

EMD Millipore Amnis ImageStreamX MkII INSPIRE User Training

Classification:


User Training

Effective Date: 01/03/2017

Revision Date: 01/25/2023

ID: UT SOP004.2

Page 3 of 4

- **Magnification:** select the 20X, 40X or 60X objective.
 - a. If viewing internal components of the cell, need to run 60X objective
 - b. 20X= 1 um/pixel
 - c. 40X= 0.5 um/pixel
 - d. 60X= 0.3 um/pixel
- **Fluidics:** Select appropriate Speed/sensitivity for your experiment
 - a. If viewing internal components of the cell, need to run on Lo speed/Hi sensitivity
- Click **Load** and place a fully-stained sample on the uptake port.
 - a. Note 15 uL is the minimum amount of liquid uptake
- Adjust the laser powers to maximize signal and prevent saturation.
 - a. Use *Raw Max Pixel_MC_ChX by Area_MChX* to determine saturation. Max output is 4000.
 - b. Be careful to closely look at axis labels as the scales will automatically update. This can be overridden by R-click Graph> Graph properties> Scaling> Manual > Maximum= 4000 (make sure scale is linear for Raw Max Pixel)
- Create plots and gates to identify the cells to collect:
 - a. *Cells in Focus: Gradient RMS_M01_Ch01*
 - b. *Singlets: Area_M01 vs. Aspect Ratio_M01*
 - c. *Phenotyping (signal intensity): use Intensity_MC_ChX*
- **Compensation** > Create Matrix... or click on the Wizards icon () and follow the prompts.
 - a. Note: Single stained samples must be collected WITHOUT brightfield/SSC channels
 - b. Collect at least 1000 positive events
 - i. Compensation wizard is looking for ONLY positive events, so it is preferable to save only events in the positive gate defined.
 - c. Run PBS between samples if worried about residual sample
 - d. Run DNA dye last
 - e. Beads are not sufficient to compensate, must use cells
- Collect all experimental samples (return the remaining sample).
- File > Save Template.

4. Shutdown

- Between users:
 - a. FACSClean (3 minutes)
 - b. diH₂O (3 minutes)
- Last user
 - a. Click **Shutdown** (sterilizes the instrument .43min)
 - b. Select “Shutdown after sterilize” (powers off all system components)
 - c. Do not exit program

Note: Sample concentration 1x10⁷-1x10⁸/mL.

Note: BF usually is set to ch1&9, DF (ch6 – SSC) should be between 3-5 mW, and single cells are visualized with a BF Area vs. Aspect Ratio.

Note: If a sample without DNA dye follows a sample with DNA dye, Load FACSClean followed by 1X PBS for a minute each.



UNC Flow Cytometry Core Facility

Title:

EMD Millipore Amnis ImageStreamX MkII INSPIRE User Training

Classification:

User Training

Effective Date: 01/03/2017

Revision Date: 01/25/2023

ID: UT SOP004.2

Page 4 of 4

5. Data Analysis Tips

- Sample Export File Types:
 - a) **.rif = raw image file**
 - b) .ctm = compensation matrix
 - c) .cif = compensated image file
 - d) **.daf = data analysis file**
 - e) .ast = tempate file
 - i) Use for batch Analysis
- F1 in IDEAS will bring up user manual; type to search

6. References

- ImageStreamX ® System Software User's Manual Version Mark II, January 2013.

7. Revisions

SOP Version Number	Date	Tracked Changes (clearly list changes made & why)	Employee
UT SOP004.2	1/25/2023	General Updates throughout to aid in user training	Ayrianna Woody