# HBE1 AIRWAY EPITHELIAL CELL LINE CULTURE RECOMMENDATIONS

- **I. SAMPLES:** Living cells at about 50% confluence in T25 flasks, filled with culture medium. The excess can be decanted, filtered to remove debris, stored at 4°C, and used over the subsequent three weeks.
- **II. Medium:** F12+7X medium (see below)
- **III. Feeding:** Change medium 2 or 3 times a week, depending on cell density. Use 5 mL for a T25 flask. Doubling time is about 36 hours.
- **IV. Passaging:** Passage when 70-90% confluent, or every 7-10 days if not that confluent. Wash cells with Phosphate Buffered Saline (PBS), then treat with 0.1% trypsin + 1 mM EDTA until cells release (usually 5-8 minutes). Neutralize trypsin with 2 fold excess Soybean Trypsin Inhibitor. Maximum expansion is 1:3, and the cells proliferate best if plated at ≥ 30% confluence. We have switched to collagen coated [*Types 1 and 3 collagen, "Purecol", Advanced Biomatrix*] provide better attachment. Proliferation is equivalent on coated and uncoated tissue culture plastic.

### V. Tips:

Passaging without expansion or with contraction is sometimes useful to stimulate growth.

Expand the lines and prepare cryopreserved stocks before starting experiments. Cryopreserve with 5% FBS and 5%DMSO, freeze at -70°C overnight, then transfer to liquid nitrogen. We freeze 5-6 vials at a time, and thaw one after three days to confirm viability.

To induce polarized differentiation, we seed permeable dishes at confluent density, and supplement medium with 50% (DME +2% FBS that has been conditioned for 3 days by confluent NIH 3T3 fibroblasts).

VI. REFERENCE: J. R. Yankaskas, et al., Am. J. Physiol. 264:c1219-c1230, 1993.

## **F-12 (7X) Media**

Explanation: F-12 (7X) media is prepared by adding 7 hormones to F-12 Media:

<b>Hormone Additive for:</b>	1 liter of F-12 Media		
T 1'	20 1		
Insulin	2.0 mL		
Transferrin (TF)	2.0 mL		
Hydrocortizone (HC)	1.0 mL		
ECGS	1.0 mL		
T3	1.0 mL		
EGF	1.0 mL		
Cholera Toxin	100 μL		

#### **Preparation of F-12 Media**

#### **MATERIALS:**

- 53.15g F-12 Nutrient Mixture (Ham) powder with L-glutamine without Sodium Bicarbonate (Gibco-21700091). Store powder at 4°C
- 50 mL 1.5 M Hepes, stored at 4°C
- 100 mL 0.714 M NaHCO<sub>3</sub>. Solution stored at 4°C
- 4 mL Gentamicin (Sigma G-1522, 50mg/mL). Solution stored at 4°C.
- 2.5mL Pen/Strep 1000X (Stock concentration 100,000 U/mL and 100 mg/mL respectively). Solution stored at -20°C

#### PROTOCOL:

- 1. Add 53.15g F-12 powder to approximately 4 liters  $H_2O$  in a 6L flask.
- 2. Stir on unheated stir plate at medium speed.
- 3. After F-12 powder has gone into solution, add 50 mL of 1.5 M Hepes and 100mL of 0.714 M NaHCO<sub>3</sub> and stir.
- 4. pH media to 7.2 using 1 M HCl or 1 M NaOH.
- 5. Add 4 mL Gentamicin and 2.5 mL Pen/Strep 1000X and stir.
- 6. Bring solution to final volume of 5 L in a graduated cylinder, with stirring, by adding deionized/distilled H<sub>2</sub>O.
- 7. Filter using a 0.2μM filter.
- 8. Store bottles at 4°C up to 4 weeks.

Note: If F12 Nutrient powder does not contain L-glutamine, add 730 mg/5 l of L-glutamine. Final concentration of L-glutamine will be 1m*M*.

# HORMONE PREPARATION AND STORAGE

Hormone/Source	Stock Prep.	Stock	Storage/	Final Conc.
		Conc.	Expires?	in F12-7X
Insulin	20 mg/4 mL dH2O	5 mg/ml	-20°C/	10 μg/mL
Sigma I-6634			3 mos.	
HC (Hydrocortisone)	1) 50 mg/13.8 mL	$10^{-2}  \text{M}$	1) -20°C/	$10^{-6}$ M
Sigma H-0888	70% ETOH		1 yr	
	2) dilute 1:10 in	$10^{-3}  \mathrm{M}$		
	PBS		2) -20C/	
			3 mos.	
ECGS (Endothelial Cell	15 mg/4 ml dH2O	3.75 mg/ml	-20°C/1	3.75 µg/ml
Growth Supplement)			mo.	
Sigma E2759				
EGF (Epidermal Growth	100 ug/4 mL dH20	25 μg/ml	-20°C/3	25 ng/ml
Factor)/ Sigma E4127			mos.	
T3 (Triiodothyronine)	1) 20 mg/5 mL	3 X 10 <sup>-5</sup> M	-20°C/3	3 X 10 <sup>-8</sup> M
Sigma T-6397	0.2M NaOH		mos.	
	2) add 5 mL dH2O			
	3) dilute 1:100 in			
	dH2O			
Transferrin	10 mg/4 mL dH2O	2.5 mg/ml	-20°C/3	5 μg/ml
Sigma T0665			mos.	
Cholera Toxin	0.5 mg/5 mL dH2O	100 μg/ml	4C°/1 yr	10 ng/ml
Sigma C-3012				