

Attune NxT Configuration

| | Mirror | Filter | Fluorochrome Examples | Detector |
|---------------------|--------|--------|------------------------------|----------|
| 561 50mW | 577LP | 585/16 | PE, dTomato, DsRed | YL1 |
| | 600LP | 620/15 | PE-Texas Red, AF594, mCherry | YL2 |
| | 650LP | 695/40 | PE-Cy5.5 | YL3 |
| | 740LP | 780/60 | PE-Cy7 | YL4 |
| 405 50mW | 417LP | 440/50 | BV421, Pacific Blue | VL1 |
| | 495LP | 512/25 | BV510 | VL2 |
| | 555LP | 603/48 | BV605 | VL3 |
| | 650LP | 710/50 | BV711 | VL4 |
| 488 50mW | | 488/10 | SSC | SSC |
| | 503LP | 530/30 | FITC, GFP, AF488 | BL1 |
| | 555LP | 590/40 | L/D Fixable Red | BL2 |
| | 650LP | 695/40 | PerCP-Cy5.5 | BL3 |
| 637 100mW | 645LP | 670/14 | APC | RL1 |
| | 690LP | 720/30 | AF700 | RL2 |
| | 740LP | 780/60 | APC-Cy7 | RL3 |