

CyAn ADP Configuration

| | Mirror | Filter | Fluorochrome Examples | Detector |
|--------------|--------|--------|-------------------------------|----------|
| 488 25mW | 95/5 | 488/10 | SSC | |
| | | 530/40 | FITC, AF488, GFP, CFSE | FL1 |
| | 545LP | 575/25 | PE, DsRed | FL2 |
| | 595LP | 613/20 | PE-TXrd, PI | FL3 |
| | 640LP | 680/30 | PE-Cy5, PerCP, 7AAD, PE-Cy5.5 | FL4 |
| | 730LP | 750LP | PE-Cy7, PE-AF750 | FL5 |
| 405 100mW | | 450/50 | PB, LD Violet, DAPI | FL6 |
| | 485LP | 560/80 | PO | FL7 |
| 643 60mW | 650LP | 665/20 | APC, AF647, To-Pro-3 | FL8 |
| | 730LP | 750LP | APC-Cy7, APC-AF750 | FL9 |

