

FACSria III

	Mirror	Filter	Fluorochrome Examples	Detector
488 100mW	655LP	695/40	PerCP-Cy5-5	A
	502LP	530/40 + 533N	FITC, GFP, AF488, CFSE, PKH67	B
		488/10	SSC	C
405 100mW	635LP	670/14	BV650	A
	595LP	610/20	BV605	B
	556LP	585/42	Pacific Orange	C
	502LP	530/30	AmCyan, CFP	D
		450/40	BV421, DAPI, Pacific Blue, L/D Violet	E
640 75mW	735LP	780/60	APC-Cy7	A
	685LP	730/45	AF700	B
		660/20	APC, AF647	C
532 120mW	735LP	780/60	PE-Cy7	A
	685LP	695/40	PE-Cy5-5	B
	635LP	670/30	7-AAD, PE-Cy5	C
	595LP	610/20	PE-Texas Red, PI, mCherry	D
		575/26	PE, DSred, dTomato	E