

Attune NxT Configuration

	Mirror	Filter	Fluorochrome Examples
561 50mW	577LP	585/16	PE, dTomato, DsRed
	600LP	620/15	PE-Texas Red, AF594, mCherry
	650LP	695/40	PE-Cy5.5
	740LP	780/60	PE-Cy7

405 50mW	417LP	440/50	BV421, Pacific Blue
	495LP	512/25	BV510
	555LP	603/48	BV605
	650LP	710/50	BV711

488 50mW		488/10	SSC
	503LP	530/30	FITC, GFP, AF488
	555LP	590/40	L/D Fixable Red
	650LP	695/40	PerCP-Cy5.5

637 100mW	654LP	670/14	APC
	690LP	720/30	AF700
	740LP	780/60	APC-Cy7

Detector
YL1
YL2
YL3
YL4

VL1
VL2
VL3
VL4

SSC
BL1
BL2
BL3

RL1
RL2
RL3