

Always filter your samples before running on any cytometer

Commercial Filters:

- [Falcon Cell Strainers](#) for 50mL tubes
- [Falcon 5mL Cell Strainer tubes](#)
- [Partec filters](#) sterile and non-sterile

Do It Yourself:

- [Blue tip and mesh filtration unit](#)
- Bulk nylon mesh filter (cut a nylon mesh to 1x1 cm pieces and placed it between a syringe and a needle) [can be sterilized with 70% ethanol]
- Tape Mesh by the yard onto a 96 well plate and pipet through

Bulk Mesh Filter:

- [Millipore](#)
- [Component Supply Company](#)
- [Nitex Screen](#)
- [Elko Filtering](#)
- [Porex](#)
- [Nylon mesh for hybridization ovens](#)

More Information:

- Purdue Cytometry List [thread 1](#)
- Purdue Cytometry List [thread 2](#)