

Attune NxT V6 Configuration

	Mirror	Filter	Fluorochrome Examples	Detector
488 50mw		488/10	SSC	SSC
	503LP	530/30	FITC, GFP, AF488	BL1
	555LP	695/40	PerCP-Cy5.5	BL2
637 100mw	654LP	670/14	APC	RL1
	690LP	720/30	AF700	RL2
	740LP	780/60	APC-Cy7	RL3
405 50 mw	417LP	450/40	Brilliant Violet 421, eFluor 450, Pacific Blue, DAPI	VL1
	495LP	525/50	eFluor 506, Brilliant Violet 510, Pacific Green	VL2
	555LP	610/20	Brilliant Violet 605, Pacific Orange	VL3
	680LP	660/20	Super Bright 645, Brilliant Violet 650	VL4
	635LP	710/50	Super Bright 702, Brilliant Violet 711	VL5
	740LP	780/60	Brilliant Violet 786	VL6
561 50 mw	577LP	585/16	PE, dTomato, DsRed	YL1
	600LP	620/15	PE-Texas Red, AF594, mCherry	YL2
	650LP	780/60	PE-Cy7	YL3