

Connect iPhone to Hearing Aids

- Better hearing for phone calls
- Better hearing for speech practice

How to connect iPhone to Hearing Aids

- 1. Open hearing aid battery doors



How to connect iPhone to Hearing Aids

- 2. Unlock iPhone
- 3. Open Settings
- 4. Tap General
- 5. Tap Accessibility
- 6. Tap Hearing Aids



How to connect iPhone to Hearing Aids

- 7. Close hearing aid battery doors




How to connect iPhone to Hearing Aids

- 8. Tap name of hearing aid
- 9. Tap "Pair"
- 10. Wait one minute



How to connect iPhone to Hearing Aids

- 11. Confirm
 - Check mark next to hearing aid
 - Control on Lock Screen



This document is based on Apple's directions to pair an iPhone to hearing aids retrieved from: <https://support.apple.com/en-us/HT201466>